



Recombinant Human Actin-related protein 5 (ACTR5)

Product Code	CSB-BP884496HU
Abbreviation	ACTR5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H9F9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAANVFPFRD ARAAPDPVLE AGPVAHGPLP VPLVLDNGSF QVRAGWACPG QDPGPEPRLQ FRAVCARGRG GARGASGPQV GNALGSLEPL RWMLRSPFDR NVPVNLELQE LLLDYSFQHL GVSSQGCVDH PIVLTEAVCN PLYSRQMMSE LLFECYGIPK VAYGIDSLFS FYHNKPKNSM CSGLIISGY QCTHVLPILE GRLDAKNCKR INLGGSQAAG YLQRLQLKY PGHLAAITLS RMEEILHEHS YIAEDYVEEL HKWRCPDYYE NNVHKMQLPF SSKLLGSTLT SEEKQERRQQ QLRRRLQELNA RRREEKLQLD QERLDRLLYV QELLEDGQMD QFHKALIELN MDSPEELQSY IQKLSIAVEQ AKQKILQAEV NLEVDVVD SK PETPDLEQLE PSLEDVESMN DFDPLFSEET PGVEKPVTTV QPVFNLAAYH QLVFGTERIR APEIIFQPSL IGEEQAGIAE TLQYILD RYP KDIQEMLVQN VFLTGGNTMY PGMKARMEKE LLEM RPFRSS FQVQLASNPV LDAWYGARDW ALNHLDDNEV WITRKEYEEK GGEYLKEHCA SNIYVPIRLP KQASRSSDAQ ASSKGSAAAGG GGAGEQA
Source	Baculovirus
Target Names	ACTR5
Protein Names	Recommended name: Actin-related protein 5 Short name= hARP5 Alternative name(s): Sarcoma antigen NY-SAR-16
Expression Region	1-607
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.