



# Recombinant Human Transcription factor COE2 (EBF2)

<b>Product Code</b>	CSB-MP884503HU
<b>Abbreviation</b>	EBF2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9HAK2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFGIQDTLGR GPTLKEKSLG AEMDSVRSWV RNVGVVDANV AAQSGVALSR AHFEKQPPSN LRKSNFFHFV LALYDRQGQP VEIERTAFVD FVENDKEQGN EKTNNGTHYK LQLLYSNGVR TEQDLVRLI DSVTKQPIAY EGQKNKPEMC RVLLTHEVMC SRCCEKKSCG NRNETPSDPV IIDRFFLKFF LKCNQNCLKT AGNPRDMRRF QVVLSTTVNV DGHVLAVSDN MFVHNNSKHG RRARRLDPSE ATPCIKAISP SEGWTTGGAM VIIIIGNFFD GLQVVFGTML VWSELITPHA IRVQTPPRHI PGVVEVTLISY KSKQFCKGAP GRFIYTALNE PTIDYGFQRL QKVIPRHPGD PERLAKEMLL KRAADLVEAL YGTPHNNQDI ILKRAADIAE ALYSVPRNPS QLPALSSSPA HSGMMGINSY GSQLGVSISE STQGNNQGYI RNTSSISPRG YSSSSTPQQS NYSTSSNSMN GYSNVPMANL GVPGSPGFLN GSPTGSPYGI MSSSPTVGSS STSSILPFSS SVFPAVKQKS AFAPVIRPQG SPSPACSSGN GNGFRAMTGL VVPPM
<b>Source</b>	Mammalian cell
<b>Target Names</b>	EBF2
<b>Protein Names</b>	Recommended name: Transcription factor COE2 Alternative name(s): Early B-cell factor 2 Short name= EBF-2
<b>Expression Region</b>	1-575
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.