



Recombinant Human Transcription factor COE2 (EBF2)

Product Code	CSB-MP884503HU
Abbreviation	EBF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HAK2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFGIQDTLGR GPTLKEKSLG AEMDSVRSWV RNVGVVDANV AAQSGVALSR AHFEKQPPSN LRKSNFFHFV LALYDRQGQP VEIERTAFVD FVENDKEQGN EKTNNGTHYK LQLLYSNGVR TEQDLYVRLI DSVTKQPIAY EGQKNKPEMC RVLLTHEVMC SRCCEKKSCG NRNETPSDPV IIDRFFLKFF LKCNQNCLKT AGNPRDMRRF QVVLSTTVNV DGHVLAVSDN MFVHNNSKHG RRARRLDPSE ATPCIKAISP SEGWTTGGAM VIIIIGDNFFD GLQVVFGTML VWSELITPHA IRVQTPPRHI PGVVEVTLISY KSKQFCKGAP GRFIYTALNE PTIDYGFQRL QKVIPRHPGD PERLAKEMLL KRAADLVEAL YGTPHNNQDI ILKRAADIAE ALYSVPRNPS QLPALSSSPA HSGMMGINSY GSQLGVSISE STQGNNQGYI RNTSSISPRG YSSSSTPQQS NYSTSSNSMN GYSNVPMANL GVPGSPGFLN GSPTGSPYGI MSSSPTVGSS STSSILPFSS SVFPAVKQKS AFAPVIRPQG SPSPACSSGN GNGFRAMTGL VVPPM
Source	Mammalian cell
Target Names	EBF2
Protein Names	Recommended name: Transcription factor COE2 Alternative name(s): Early B-cell factor 2 Short name= EBF-2
Expression Region	1-575
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.