



Recombinant Human Vasculin-like protein 1 (GPBP1L1)

Product Code	CSB-BP884517HU
Abbreviation	GPBP1L1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9HC44
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAQHDFVPAW LNFSTPQSAK SPTATFEKHG EHLPRGEGRF GVSRRRHNSG DGFFNNGPLR TAGDSWHQPS LFRHDSVDSG VSKGAYAGIT GNPSGWHSSS RGHGMSQRS GGGTGNHRHW NGSFHSRKGC AFQEKPPMEI REEKEDKVE KLQFEEEDFP SLNPEAGKQH QPCRPIGTPS GVWENPPSAK QPSKMLVIKK VSKEDPAAAF SAAFTSPGSH HANGNKLSSV VPSVYKNLVP KPVPPPSKPN AWKANRMEHK SGSLSSSRES AFTSPISVTK PVVLASGAAL SSPKESPSST TPPIEISSSR LTKLTRRTTD RKSEFLKTLK DDRNGDFSEN RDCDKLEDLE DNSTPEPKEN GEEGCHQNGL ALPVVEEGEV LSHSLEAEHR LLKAMGWQEY PENDENCLPL TEDELKEFHM KTEQLRRNGF GKNGFLQSRS SSLFSPWRST CKAEFEDSDT ETSSSETSDD DAWK
Source	Baculovirus
Target Names	GPBP1L1
Protein Names	Recommended name: Vasculin-like protein 1 Alternative name(s): GC-rich promoter-binding protein 1-like 1
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.