



Recombinant Rat Nucleobindin-2 (Nucb2)

Product Code	CSB-YP884912RA
Abbreviation	Nucb2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9JI85
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	VPIDVD KTKVHNVEPV ESARIEPPDT GLYYDEYLKQ VIEVLETDPH FREKLQKADI EEIRSGRLSQ ELDLVSHKVR TRLDELKRQE VGRLRMLIKA KLDALQDTGM NHHLLLKQFE HLNHQNPDTF ESKDLMLIK AATADLEQYD RTRHEEFKKY EMMKEHERRE YLKTLSEEKR KEEEAKFAEM KRKHEDHPKV NHPGSKDQLK EVWEETDGLD PNDFDPKTF KLHDVNNDGF LDEQEALF TKELDKVYNP QNAEDDMIEM EEERLRMREH VMNEIDNNKD RLVTL EEFLR ATEKKEFLEP DSWETLDQQQ LFTEEELKEY ESIIAIQESE LKKKADELQK QKEELQRQHD HLEAQKQEYQ QAVQQLEQKK FQQGIAPSGP AGELKFEPHT
Source	Yeast
Target Names	Nucb2
Protein Names	Recommended name: Nucleobindin-2 Alternative name(s): DNA-binding protein NEFA
Expression Region	25-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>