



Recombinant Mouse Mothers against decapentaplegic homolog 9 (Smad9)

Product Code	CSB-MP884919MO
Abbreviation	Smad9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9JIW5
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MHPSTPISSL FSFTSPAVKR LLGWKQGDEE EKWAEKAVDS LVKKLKKKKG AMDELERALS CPGQPSKCVT IPRSLDGRLQ VSHRKGLPHV IYCRVWRWPD LQSHHELKPL ECCEFPFGSK QKEVCINPYH YRRVETPVLP PVLVPRHSEY NPQLSLLAKF RSASLHSEPL MPH NATYPDS FQQSLCPAPP SSPGHVFPQS PCPTSYPHSP GSPSESDSPY QHSDFRPVCY EEPQHWCSVA YYELNNRVGE TFQASSRSVL IDGFTDPSNN RNRFLGLLS NVNRNSTIEN TRRHIGKGVH LYYVGGEVYA ECVSDSSIFV QSRNCNYQH G FHPATVCKIP SGCSLKVFNN QLFAQLLAQS VHHGFEVVEY LTKMCTIRMS FVKGWGAEYH RQDVTSTPCW IEIHLHGPLQ WLDKVL TQMG SPHNPISSVS
Source	Mammalian cell
Target Names	Smad9
Protein Names	Recommended name: Mothers against decapentaplegic homolog 9 Short name= MAD homolog 9 Short name= Mothers against DPP homolog 9 Alternative name(s): SMAD family member 9 Short name= SMAD 9 Short name= Smad9 Smad8
Expression Region	1-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.