



Recombinant Arabidopsis thaliana Phosphoglycerate kinase 1, chloroplastic (PGK1)

Product Code	CSB-MP885315DOA
Abbreviation	PGK1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LD57
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	SMAKK SVGDLTSADL KGKKVFVRAD LNVPLDDNQT ITDDTRIRAA IPTIKYLIEN GAKVILSTHL GRPKGVTPKF SLAPLVPRLS ELLGIEVTKA DDCIGPEVES LVASLPEGGV LLENVRFYK EEEKNDPEFA KKLASLADLY VNDAFGTAHR AHASTEGVTK FLKPSVAGFL LQKELDYLVG AVSNPKRPFA AIVGGSKVSS KIGVIESLLE KCDILLGGG MIFTFYKAQG LSVGSSLVEE DKLELATELL AKAKAKGVSL LLPTDVVVAD KFAPDANSKI VPASGIEDGW MGLDIGPDSI KTFNEALDTT QTVIWNGPMG VFEMEKFAAG TEAIANKLAE LSEKGVTTII GGGDSVAAVE KVG VAGVMSH ISTGGGASLE LLEGKVLPGV IALDEAIPVT V
Source	Mammalian cell
Target Names	PGK1
Protein Names	Recommended name: Phosphoglycerate kinase 1, chloroplastic EC= 2.7.2.3
Expression Region	76-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.