



Recombinant Hevea brasiliensis Enolase 1 (ENO1)

Product Code	CSB-MP885326HWI
Abbreviation	ENO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LEJ0
Product Type	Recombinant Protein
Immunogen Species	Hevea brasiliensis (Para rubber tree) (Siphonia brasiliensis)
Purity	≥85% (SDS-PAGE)
Sequence	MAITIVSVRA RQIFDSRGNP TVEADV KLS D GYLARAAVPS GASTGIYEAL ELRDGGSDYL GKGVS KAVEN VNIIGPALV GKDPTDQVGI DNFMVQQLDG TVNEWGWCKQ KLGANAILAV SLAVCKAGAH VKGIPLYEHI ANLAGNKNLV LPVPAFNVIN GGSHAGNKLA MQEFMILPVG ASSFKEAMKM GAEVYHHLKS VIKKKYGQDA TNVGDEGGFA PNIQENKEGL ELLKTAIAKA GYTGVVIGM DVAASEFYGS DQTYDLNFKE ENNNGSQKIS GEALKDLYKS FVAEYPIVSI EDPFDQDDWA HYAKLTSEIG EKQIVGDDL LVTNPKRVEK AIKEKACNAL LLKVNQIGSV TESIEAVKMS KRAGWGVMS HRSGETEDTF IADLSVGLAT GQIKTGAPCR SERLAKYNQL LRIEEELGSE AVYAGANFRK PVEPY
Source	Mammalian cell
Target Names	ENO1
Protein Names	Recommended name: Enolase 1 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 1 2-phosphoglycerate dehydratase 1 Allergen= Hev b 9
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.