



Recombinant Arabidopsis thaliana Ran-binding protein 1 homolog a (RANBP1A)

Product Code	CSB-MP885378DOA
Abbreviation	RANBP1A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LMK7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MATNEPEHEH RDEEEAGANE DEDTGAQVAP IVRLEEVAVT TGEEDEDAVAL DLKSKLYRFD KDANQWKERG AGTVKFLKHK NTGKIRLVMR QSKTLKICAN HFVKSGMSVQ EHVGNKSCV WHARDFADGE LKDELFCIRF ASIENCKTFM QKFKEVAESE EEKEESKDA DTAGLLEKLT VEETKTEEKT EAKAVETAKT EVKAEKKES EAEKSGEAKK TEESGPST
Source	Mammalian cell
Target Names	RANBP1A
Protein Names	Recommended name: Ran-binding protein 1 homolog a Alternative name(s): Ran-binding protein siRanBP
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.