



Recombinant Arabidopsis thaliana F-box protein At1g48060 (At1g48060)

Product Code	CSB-MP885385DOA
Abbreviation	At1g48060
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LNF6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKPQEEEEKN ENMARKRSKS SSSLSIPLDI ATDIFLRLPA KSVVRFSCVA KHWSSITTAP YFTNSFETRP NLLFFFKEAD RFFVVTVPKP NRRPNESVSH TSSQILDSYQ TPYPKHSCFT IKTESVHGLI CFQRGTKPIV WNPTMRKFKP LRKPDKSWES LTVFLGYDPL ERTHKVVAMP CDKASDECRI LTLGSADQES WRTVKNTYKH HPSRGNRKNN YGPCRCINGV LYYLAEIDRH RLIGKFAFLG HSNVKNSLPM WVLEGEKKGE WSTYNFLPLS HYDRRLEFHF KLIGITNDGE LIYVPNTVFK TFEVIYIDPI RKTFSLVKYE GVADKGFQRQR NGLGEDKPLR GIQYSPNHVE TLLSL
Source	Mammalian cell
Target Names	At1g48060
Protein Names	Recommended name: F-box protein At1g48060
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.