



Recombinant Arabidopsis thaliana Peroxisomal (S)-2-hydroxy-acid oxidase GLO2 (GLO2)

Product Code	CSB-MP885413DOA
Abbreviation	GLO2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LRS0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEITNVTEYD AIAKAKLPKM VYDYYASGAE DQWTLQENRN AFARILFRPR ILIDV NKIDM ATTVLGFKIS MPIMVAPTAF QKMAHPDGEY ATARAASAAG TIMTLSSWAT SSVEEVASTG PGIRFFQLYV YKNRKVVEQL VRRAEKAGFK AIALTVDTPR LGRRES DIKN RFTLPPNLT LKNFEGLDLGK MDEANDSGLA SYVAGQIDRT LSWKDIQWLQ TITNMPILVK GVL TGEDARI AIQAGAAGII VSNHGARQLD YVPATISALE EVVKATQGRV PVFLDGGVRR GTDVF KALAL GASGIFIGRP VVFALAAEGE AGVKKVLQML RDEFELTMAL SGC RSLSEIT RNHIVTEWDT PRHLPRL
Source	Mammalian cell
Target Names	GLO2
Protein Names	Recommended name: Peroxisomal (S)-2-hydroxy-acid oxidase GLO2 EC= 1.1.3.15 Alternative name(s): Glycolate oxidase 1 Short name= AtGLO2 Short name= GOX 1 Short chain alpha-hydroxy acid oxidase GLO2
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.