



Recombinant *Oryza sativa* subsp. japonica Proteasome subunit alpha type-3 (PAG1)

Product Code	CSB-EP885426OFG
Abbreviation	PAG1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LSU0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSSIGTGYDL SVTTFSPDGR VFQVEYATKA VDNSGTVVGI KCKDGIVLGV EKLVTSKMML EGSNRRIVSV HWHSGGLAVAG LAADGRQIVS RAKSEAASYE KVYGEPISVK ELADRVASYV HLCTLYWWLR PFGCGVILGG YDRDGPQLYM IEPSGVSYKY FGAALGKGRQ AAKTEIEKLG LSELTCREGI VEVAKIIYGV HDEAKDKAFE LELSWICDES NRQHQKVPAD LLEQAKVAAQ AALEEMDAD
Source	E.coli
Target Names	PAG1
Protein Names	Recommended name: Proteasome subunit alpha type-3 EC= 3.4.25.1 Alternative name(s): 20S proteasome alpha subunit G 20S proteasome subunit alpha-7
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.