



# Recombinant Arabidopsis thaliana Ethylene-responsive transcription factor RAP2-2 (RAP2-2)

<b>Product Code</b>	CSB-BP885447DOA
<b>Abbreviation</b>	RAP2-2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LUM4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MCGGAIISDF IPPRSLRVT NEFIWPD LKN KVKASKKRSN KRSDFFDLDD DFEADFQGFK DDSAFDCEDD DDVFNVKPF VFTATTKPVA SAFVSTGIYL VGSAYAKKTV ESAEQAEKSS KRKRKNQYRG IRQRPWGKWA AEIRDPRKGS REWLGTFDTA EEAARAYDAA ARRIRGTKAK VNFPEEKNSP VVSQKRPSAK TNNLQKSVAK PNKSVTLVQQ PTHLSQQYCN NSF DNSFGDM SFMEEKPMY NNQFGLTNSF DAGGNGYQY FSSDQGSNSF DCSEFGWSDH GPKTPEISSM LVNNNEASV EETNAAKLK PNSDESDDL M AYLDNALWDT PLEVEAMLGA DAGAVTQEEE NPVELWSLDE INFMLEGDF
<b>Source</b>	Baculovirus
<b>Target Names</b>	RAP2-2
<b>Protein Names</b>	Recommended name: Ethylene-responsive transcription factor RAP2-2 Alternative name(s): Protein RELATED TO APETALA2 2
<b>Expression Region</b>	1-379
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.