



Recombinant Arabidopsis thaliana F-box/kelch-repeat protein At3g16580 (At3g16580)

Product Code	CSB-MP885450DOA
Abbreviation	At3g16580
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LUS6
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAHEEKRPE FLSLSPWELI EEILSRVPE SLLRFKTVSK QWNALFRDKT FINNHKMTFR FILATKSKIY SVSIDPKIVV RELTLDIPGL ESHEIPKKLV DCDKLLLCDM EKGVVLWNPW LRHSTWIDQG SNHTRMESYG IGYNNKGSYK IFAFCDRKEN HTQRLLTIHD SASDAWKDRE PIDNSQGKQI VHNIYTKISG VSLNGNLYLV TYFETDLDVY HLIENSSSE SVVKFCDLPC GTSNFLKDAF VLRVFEGDRF SLLKQCHATK KIEIWSKYK INNNLDRDVE WIKFMEVSSP NLPDLVDGFD SQPSYFIEDK RLVVCSCNET GRAWIYVFGE NKLISKTQID SVVDLWPSHW TFIPSLVPVP RAQREPEAEL QV
Source	Mammalian cell
Target Names	At3g16580
Protein Names	Recommended name: F-box/kelch-repeat protein At3g16580
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.