



Recombinant Arabidopsis thaliana B3 domain-containing protein At5g60130 (At5g60130)

Product Code	CSB-MP885457DOA
Abbreviation	At5g60130
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LVG1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNKGVSLDNC VPKFFKVYLP DESGDDMVLP ISFTRFLPKS LPETVTVRSI SGNIWKLELK KCCGDDDETEK FVMVNGWKRI VKDEDLKGGD FLEFEFDGSS CFHFCIYEHR TMCKKIRRSS DQSEEIKVES DSDEQNQASD DVLSLDEDDD DSDYNCGEDN DSDDYADEAA VEKDDNDADD EDVDNVADDV PVEDDDYVEA FDSRDHAKAD DDEDERQYL DDRENPSFTL ILNPKKKSQ LIPARVIKDY DLHFPESITL VDPLVKKFGT LEKQIKIQTN GSVFVKGFSG IIRRNKVKTT DKMIFEIKKT GDNNLVQTIK IHIIISG
Source	Mammalian cell
Target Names	At5g60130
Protein Names	Recommended name: B3 domain-containing protein At5g60130
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.