



Recombinant Arabidopsis thaliana F-box/LRR-repeat protein At3g59250 (At3g59250)

Product Code	CSB-BP885469DOA
Abbreviation	At3g59250
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LX46
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MDIGPKDKIS NLPEALICHI LSFLPIEDSA LTSVLSKRWR YLFAFRPNLV FDDSVCLRPP CGHGVIPDYV TSWILNVLER GVLDDLQVN VNGVRLPSKV FVSKSLVRLR IESGNVRGID VDDVFLPKL SLYLHTIMLG KGEDCFEKL SGCHVLEELV LNNVYSNVWN RSVSSKTLKR LTLWCTESYD KNPHSVLFDT PNLVYLKYSYD YVAAKYPKVN FSSLVDATIG LAMTSDQDAN ASYDALVGNA TDFLMGICNV QILCLSANSL AVLTFCCESI PVFNLIQLT IKTNQSVGWE SLPALLKNCP ILETLVFEGL LHMYTIKCGD VDGCLCKYSG EVPSCLSSSP VKVLKILRFG EYGIVNQIEM IKIFLETMPH LEQLTVYYNT PIGGDLIEEV SSQLKNFARE ASPNCEIQVI PDVLS
Source	Baculovirus
Target Names	At3g59250
Protein Names	Recommended name: F-box/LRR-repeat protein At3g59250
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.