



# Recombinant Arabidopsis thaliana F-box/LRR-repeat protein At3g59250 (At3g59250)

<b>Product Code</b>	CSB-EP885469DOA-B
<b>Abbreviation</b>	At3g59250
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9LX46
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDIGPKDKIS NLPEALICHI LSFLPIEDSA LTSVLSKRWR YLFAFRPNLV FDDSVCLRPP CGHGVIPDYV TSWILNVLER GVLDDLQVN VNGVRLPSKV FVSKSLVRLR IESGNVRGID VDDVFLPKL SLYLHTIMLG KGEDCFEKL SGCHVLEELV LNNVYSNVWN RSVSSKTLKR LTLWCTESYD KNPHSVLFDT PNLVYLYKYS YVAAKYPKVN FSSLVDATIG LAMTSDQDAN ASYDALVGNA TDFLMGICNV QILCLSANSL AVLTFCCESI PVFNLIQLT IKTNQSVGWE SLPALLKNCP ILETLVFEGL LHMYTIKCGD VDGCLCKYSG EVPSCLSSSP VKVLKILRFG EYGIVNQIEM IKIFLETMPH LEQLTVYYNT PIGGDLIEEV SSQLKNFARE ASPNCEIQVI PDVLS
<b>Source</b>	E.coli
<b>Target Names</b>	At3g59250
<b>Protein Names</b>	Recommended name: F-box/LRR-repeat protein At3g59250
<b>Expression Region</b>	1-425
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.