



Recombinant Arabidopsis thaliana Putative F-box/LRR-repeat protein At3g58880 (At3g58880)

Product Code	CSB-BP885476DOA
Abbreviation	At3g58880
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9LXR4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVDLVSSLPD DLLGHILSLL TTKEAALTSI LSKRWRYLIA FVPYLEFDDS AFLNPEEGKQ TREGTRQSFI DFVDRVLALH GDSPIRKFSL KCKTGVDL DL LNQWICNVLQ RGVLLIDLSM DLGHRMCFIE IFMSRTLVEL KLGSGCRIA F GPEHISALPM LKTLTLDSVS WSDSGQLERL LSACPAL EAL NLANVHGSYP NATVSIASLK TLTIKSVSLS GPAHVFSFD TPNLLCLNYTA LFEDDYPLVN LEYLVEAQIK FVLTDRLIKL VSVRKNGLLM LSEVQKLIRG ISSVRKLYLS PGTLQVLGQC SQAMPVFNNL TFLVISSMD IRWQAMPVLL KNCPRLETLV IKGGLVHCVA ADCGDACTCI SREEKGRSLA SCPVKRLEIR EFQGTREME MIKHFYYCFL CLKEMEIYVK DGRPQFLAPL ASDSNPFKHL LGRNVIKVVH GSLT
Source	Baculovirus
Target Names	At3g58880
Protein Names	Recommended name: Putative F-box/LRR-repeat protein At3g58880
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.