



Recombinant Arabidopsis thaliana Dof zinc finger protein DOF4.1 (DOF4.1)

Product Code	CSB-BP885510DOA
Abbreviation	DOF4.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M161
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDHHQYHHHD QYQHQMSTST NNNSYNTIVT TQPPPTTTM DSTTATTMIM DDEKLMSTM STRPQEPRNC PRCNSSNTKF CYYNNYSLAQ PRYLCKSCRR YWTEGGSLRN VPVGGGSRKN KKLFPNSST SSSTKNLPDL NPPFVFTSSA SSSNPSKTHQ NNNDLSLSFS SPMQDKRAQG HYGHFSEQVV TGGQNCLFQA PMGMIQFRQE YDHEHPKKNL GFSLDRNEEE IGNHDFVFN EEGSKMMYPY GDHEDRQQHH HVRHDDGNKK REGGSSNELW SGIILGGDSG GPTW
Source	Baculovirus
Target Names	DOF4.1
Protein Names	Recommended name: Dof zinc finger protein DOF4.1 Short name= AtDOF4.1
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.