



Recombinant Arabidopsis thaliana Branched-chain-amino-acid aminotransferase 2, chloroplastic (BCAT2)

Product Code	CSB-MP885535DOA
Abbreviation	BCAT2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9M439
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	LQTFAYKN AQAASALREE RKKPLYQNGD DVYADLDWDN LGFGLNPADY MYVMKCSKDG EFTQGELSPY GNIQLSPSAG VLNYGQAIYE GTKAYRKENG KLLLFRPDHN AIRMKLGAER MLMPSPSVDQ FVNAVKQTAL ANKRWWPPAG KGTLYIRPLL MGSGPILGLG PAPEYTFIVY ASPVGNFYFKE GMAALNLYVE EEYVRAAPGG AGGVKSITNY APVLKALSRA KSRGFSDVLY LDSVKKKYLE EASSCNV FVV KGRTISTPAT NGTILEGITR KSVMEIASDQ GYQVVEKAVH VDEVMDADEV FCTGTAVVVA PVGTITYQEK RVEYKTGDES VCQKLRSLV GIQTGLIEDN KGWVTDIN
Source	Mammalian cell
Target Names	BCAT2
Protein Names	Recommended name: Branched-chain-amino-acid aminotransferase 2, chloroplastic Short name= Atbcac-2 EC= 2.6.1.42
Expression Region	23-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.