



# Recombinant Arabidopsis thaliana 26S protease regulatory subunit S10B homolog B (RPT4B)

<b>Product Code</b>	CSB-MP885565DOA
<b>Abbreviation</b>	RPT4B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9MAK9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSDGDDAARR RTAAVTDYRK KLLHHKELES RVRTARENLR AAKKEFNKTE DDLKSLQSVG QIIGEVLRPL DNERLIVKAS SGPRYVVGCR SKVDKEKLTS GTRVVLDMTT LTIMRALPRE VDPVVYNMLH EDPGNISYSA VGGLGDQIRE LRESIELPLM NPELFLRVGI KPPKGVLLYG PPGTGKTLA RAIASNIDAN FLKVSSAII DKYIGESARL IREMFNYARE HQPCIIFMDE IDAIGGRRFS EGTSADREIQ RTLMELLNQL DGFDQLGKVK MIMATNRPDV LDPALLRPGR LDRKIEIPLP NEQSRMEILK IHASGIAKHG EIDYEAIVKL GEGFNGADLR NICTEAGMFA IRAERDYVIH EDFMKAVRKL SEAKKLESSS HYNADFGKE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RPT4B
<b>Protein Names</b>	Recommended name: 26S protease regulatory subunit S10B homolog B Alternative name(s): 26S proteasome AAA-ATPase subunit RPT4b 26S proteasome subunit S10B homolog B Regulatory particle triple-A ATPase subunit 4b
<b>Expression Region</b>	1-399
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.