



# Recombinant Rabbit Filamin-B (FLNB), partial

<b>Product Code</b>	CSB-BP885623RB
<b>Abbreviation</b>	FLNB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9MZD2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GTRLNGAKGK IDAKVHSPSG AVEECHVSEL EPDKYAVRFI PHENGIHTID VKFNGSHVVG SPFKVRVGEF GQAGNPALVS AYGAGLEGGT TGIQSEFFIN TTRAGPGTLS VTIEGPSKVK MDCQETPEGY KVMYTPMAPG NYLIGVKYGG PNHIVGSPFK AKVTGQRLVG PGSTNETSSI LVESVTRSST ETCYSAIPKA SSDASKVTSK GAGLSKAFVG QKSSFLVDCS KAGSNMLLIG VHGPPTPCEE VSMKHVGSQQ YNVTYVVKER GEYVLAVKWG EEHIPGSPFH VTVP
<b>Source</b>	Baculovirus
<b>Target Names</b>	FLNB
<b>Protein Names</b>	Recommended name: Filamin-B Short name= FLN-B Alternative name(s): ABP-280-like protein Actin-binding-like protein Beta-filamin C254
<b>Expression Region</b>	1-294
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Target Details</b>	This gene encodes a member of the filamin family. The encoded protein interacts with glycoprotein Ib alpha as part of the process to repair vascular injuries. The platelet glycoprotein Ib complex includes glycoprotein Ib alpha, and it binds the actin cytoskeleton. Mutations in this gene have been found in several conditions: atelosteogenesis type 1 and type 3; boomerang dysplasia; autosomal dominant Larsen syndrome; and spondylocarpotarsal synostosis syndrome. Multiple alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## **Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.