



Recombinant Pan paniscus B-cell receptor CD22 (CD22)

Product Code	CSB-EP885632EQU
Abbreviation	CD22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9N1E5
Product Type	Recombinant Protein
Immunogen Species	Pan paniscus (Pygmy chimpanzee) (Bonobo)
Purity	>85% (SDS-PAGE)
Sequence	D SSKWAFEHPE TLYAWEGACV WIPCTYRALD RDLESFILFH NPEYNKNTSK FDGTRLYEST EDGKVPSEQK RVQFLGDKNK NCTLSIHPVH VNDSGQLGLR MESKTAKWME RIHLNVSERP FPPHIQLPPE IQESQEVTLT CLLNFSCYGY PIQLQWLLEG VPMRQAAVTS TSLTIKSVFT RSEKLFSPQW SHHGKIVTCQ LQDADGKFLS NDTVQLNVKH TPKLEIKVTP SDAIVREGES VTMTCEVSSS NPEYTTISWL KDGTSLLKKQN TLMLNLHEVT KDQSGKYCCQ VSNDVGPGRS AEVFLQVQYA PE
Source	E.coli
Target Names	CD22
Protein Names	Recommended name: B-cell receptor CD22 Alternative name(s): Sialic acid-binding Ig-like lectin 2 Short name= Siglec-2 CD_antigen= CD22
Expression Region	20-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.