



Recombinant Human Doublesex- and mab-3-related transcription factor 3 (DMRT3)

Product Code	CSB-EP885692HU-B
Abbreviation	DMRT3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NQL9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNGYGSPYLY MGGPVSQPPR APLQRTPKCA RCRNHGVLWSW LKGHKRYCRF KDCTCEKCIL IIERQRVMAA QVALRRQQAN ESLES LIPDS LRALPGPPPP GDAVAAPQPP PASQPSQPQP PRPAAELAAA AALRWTAEPQ PGALQAQLAK PDLTEERLGD GKSADNTEVF SDKDTDQRSS PDVAKSKGCF TPESPEIVSV EEGGYAVQKN GGNPESRPDS PKCHAEQNH LIEGPSGTVS LPFSLKANRP PLEV LKKIFP NQKPTVLELI LKGC GGD LVS AVEVLLSSRS SVTGAERTSA EPESLALPSN GHIFEHTLSS YPISSSKWSV GSAFRVPDTL RFSADSSNVV PSPLAGPLQP PFPQPPRYPL MLRNTLARSQ SSPFLPNDVT LWNTMTLQQQ YQLRSQYVSP FPSNSTSVFR SSPVLPARAT EDPRISIPDD GCPFVSKQSI YTEDDYDERS DSSDSRTLNT SS
Source	E.coli
Target Names	DMRT3
Protein Names	Recommended name: Doublesex- and mab-3-related transcription factor 3
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.