



Recombinant Human Endothelial zinc finger protein induced by tumor necrosis factor alpha (ZNF71)

Product Code	CSB-EP885697HU
Abbreviation	ZNF71
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NQZ8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKELDPKNDI SEDKLSVVG E ATGGPTRNGA RGP GSEGVWE PGSWPERPRG DAGAEWEPLG IPQGNKLLGG SVPACHELKA FANQGCVLVP PRLDDPTEKG ACPPVRRGKN FSSTSDL SKP PMPCEEK KTY DCSECGKAFS RSSSLIKHQR IHTGKPFEC DTCGKH FIER SSLTIHQ RVH TGEKPYACGD CGKAFSQRMN LTVHQRTHTG EKP YVCDVCG KAFRKTSSLT QHERIHTG EK PYACGDCGKA FSQNMHLIVH QRTHTG EKPY VCPECGRAFS QNMHLTEHQR THTG EKPYAC KECGKAFNKS SSLTLHQ RNH TGEKPYVCGE CGKAFSQSSY LIHQRFHIG VKPFECSECG KAFSKNSSLT QHQRIHTG EK PYECYICKKH FTGRSSLIVH QIVHTG EKPY VCGECGKAFS QSAYLIEHQR IHTG EKPYRC GQCGKSF IKN SSLTVHQRIH TGEKPYRCGE CGKTF SRNTN LTRHLRIHT
Source	E.coli
Target Names	ZNF71
Protein Names	Recommended name: Endothelial zinc finger protein induced by tumor necrosis factor alpha Alternative name(s): Zinc finger protein 71
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.