



# Recombinant Human Endothelial zinc finger protein induced by tumor necrosis factor alpha (ZNF71)

<b>Product Code</b>	CSB-MP885697HU
<b>Abbreviation</b>	ZNF71
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NQZ8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKELDPKNDI SEDKLSVVG E ATGGPTRNGA RGP GSEGVWE PGSWPERPRG DAGAEWEPLG IPQGNKLLGG SVPACHELKA FANQGCVLVP PRLDDPTEKG ACPPVRRGKN FSSTSDLSKP PMPCEEKKT Y DCSECGKAFS RSSSLIKHQR IHTGEKPFEC DTCGKHFIER SSLTIHQRVH TGEKPYACGD CGKAFSQRMN LTVHQRTHTG EKPYVCDVCG KAFRKTSSLT QHERIHTGEK PYACGDCGKA FSQNMHLIVH QRTHTGEKPY VCPECGRAFS QNMHLTEHQR THTGEKPYAC KECGKAFNKS SSLTLHQRNH TGEKPYVCGE CGKAFSQSSY LIHQRFHIG VKPFECSECG KAFSKNSSLT QHQRIHTGEK PYECYICKKH FTGRSSLIVH QIVHTGEKPY VCGECGKAFS QSAYLIEHQR IHTGEKPYRC GQCGKSF IKN SSLTVHQRIH TGEKPYRCGE CGKTFSRNTN LTRHLRIHT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF71
<b>Protein Names</b>	Recommended name: Endothelial zinc finger protein induced by tumor necrosis factor alpha Alternative name(s): Zinc finger protein 71
<b>Expression Region</b>	1-489
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.