



# Recombinant Human N-lysine methyltransferase SMYD2 (SMYD2)

<b>Product Code</b>	CSB-EP885705HU
<b>Abbreviation</b>	SMYD2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NRG4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRAEGLGGLE RFCSPGKGRG LRALQPFQVG DLLFSCPAYA YVLTVNERGN HCEYCFTRKE GLSKCGRCKQ AFYCNVECQK EDWPMHKLEC SPMVVFGENW NPSETVRLTA RILAKQKIHP ERTPSEKLLA VKEFESHLDK LDNEKKDLIQ SDIAALHHFY SKHLGFPDND SLVVLFAQVN CNGFTIEDEE LSHLGSAIFP DVALMNHSCC PNVIVTYKGT LAEVRAVQEI KPGEEVFVTSY IDLLYPTEDR NDRLRDSYFF TCECQECTTK DKDKAKVEIR KLSDPPKAEA IRDMVRYARN VIEEFRAKH YKSPSELLEI CELSQEKMSS VFEDSNVYML HMMYQAMGVC LYMQDWEGAL QYGQKIIKPY SKHYPLYSLN VASMWLKLGR LYMGLEHCAA GEKALKKAAIA IMEVAHGKDH PYISEIKQEI ESH
<b>Source</b>	E.coli
<b>Target Names</b>	SMYD2
<b>Protein Names</b>	Recommended name: N-lysine methyltransferase SMYD2 EC= 2.1.1.- Alternative name(s): HSKM-B Histone methyltransferase SMYD2 EC= 2.1.1.43 Lysine N-methyltransferase 3C SET and MYND domain-containing protein 2
<b>Expression Region</b>	1-433
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.