



Recombinant Human Olfactomedin-like protein 3 (OLFML3)

Product Code	CSB-YP885708HU
Abbreviation	OLFML3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NRN5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	QQHHLVEYM ERRLAALER LAQCQDQSSR HAAELRDFKN KMLPLLEVAE KEREALRTEA DTISGRVDRL EREVDYLETQ NPALPCVEFD EKVTGGPGTK GKGRRNEKYD MVTDCGYTIS QVRSMKILKR FGGPAGLWTK DPLGQTEKIY VLDGTQNDTA FVFPRLRDFT LAMAARKASR VRVFPWPVGT GQLVYGGFLY FARRPPGRP GGGEMENTLQ LIKFHLANRT VVDSSVFP AE GLIPPYGLTA DTYIDLAADE EGLWAVYATR EDDRHLCLAK LDPQTLDEQ QWDTPCPREN AEAAFVICGT LYVVYNTRPA SRARIQCSFD ASGTLTPERA ALPYFPRRYG AHASLRYNPR ERQLYAWDDG YQIVYKLEMR KKEEEV
Source	Yeast
Target Names	OLFML3
Protein Names	Recommended name: Olfactomedin-like protein 3 Alternative name(s): HNOEL-iso hOLF44
Expression Region	22-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.