



Recombinant Human SAM domain-containing protein SAMSN-1 (SAMSN1)

| | |
|--------------------------|---|
| Product Code | CSB-EP885722HU |
| Abbreviation | SAMSN1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9NSI8 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MLKRKPSNVS EKEKHQKPKR SSSFGNFDRF RNNSLSKPDD STEAHEGDPT NGSGEQSKTS NNGGGLGKKM RAISWTMKKK VGKKYIKALS EEKDEEDGEN AHPYRNSDPV IGTHTKEVSL KASDSMDSLY SGQSSSSGIT SCSDGTSNRD SFRLDDDPY SGPFCGRARV HTDFTSPYD TDSLKIKKGD IIDIICKTPM GMWTGMLNNK VGNFKFIYVD VISEEEAAPK KIKANRRSNS KKSRTLQEFL ERIHLEQYTS TLLNGYETL EDLKDIKESH LIELNIENPD DRRRLLSAAE NFLEEEIIQE QENEPEPLSL SSDISLNKSQ LDDCPRDSGC YISSGNSDNG KEDLESENLS DMVHKIITE PSD |
| Source | E.coli |
| Target Names | SAMSN1 |
| Protein Names | Recommended name: SAM domain-containing protein SAMSN-1 Alternative name(s): Hematopoietic adaptor containing SH3 and SAM domains 1 Nash1 SAM domain, SH3 domain and nuclear localization signals protein 1 SH3-SAM adaptor protein |
| Expression Region | 1-373 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.