



Recombinant Human SAM domain-containing protein SAMSN-1 (SAMSN1)

Product Code	CSB-YP885722HU
Abbreviation	SAMSN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NSI8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLKRKPSNVS EKEKHQKPKR SSSFGNFDRF RNNSLSKPDD STEAHEGDPT NGSGEQSKTS NNGGGLGKKM RAISWTMKKK VGKKYIKALS EEKDEEDGEN AHPYRNSDPV IGTHTKVS L KASDSMDSLY SGQSSSSGIT SCSDGTSNRD SFRLDDDGYPY SGPFCGRARV HTDFTSPSYD TDSLKIKKGD IIDIICKTPM GMWTGMLNNK VGNFKFIYVD VISEEEAAPK KIKANRRSNS KKSRTLQEFL ERIHLEQYTS TLLLNGYETL EDLKDIKESH LIELNIENPD DRRRLLSAAE NFLEEEIIQE QENEPEPLSL SSDISLNKSQ LDDCPRDSCG YISSGNSDNG KEDLESENLS DMVHKIITE PSD
Source	Yeast
Target Names	SAMSN1
Protein Names	Recommended name: SAM domain-containing protein SAMSN-1 Alternative name(s): Hematopoietic adaptor containing SH3 and SAM domains 1 Nash1 SAM domain, SH3 domain and nuclear localization signals protein 1 SH3-SAM adaptor protein
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.