



Recombinant Human Ufm1-specific protease 2 (UFSP2)

Product Code	CSB-YP885731HU
Abbreviation	UFSP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NUQ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVISESMDIL FRIRGGLDLA FQLATPNEIF LKKALKHVLS DLSTKLSSNA LVFRICHSSV YIWPSSDINT IPGELTDASA CKNILRFIQF EPEEDIKRKF MRKKDKKLSLSD MHQIVNIDLM LEMSTSLAAV TPIIERESGG HHYVNMTLPV DAVISVAPEE TWGKVRKLLV DAIHNQLTDM EKCILKYMKG TSIVVPEPLH FLLPGKKNLV TISYPSGIPD GQLQAYRKEL HDLFLNPHDR PYFKRSNAYH FPDEPYKDG Y IRNPHTYLNP PNMETGMIYV VQGIYGYHHY MQDRIDDNGW GCAYRSLQTI CSWFKHQGYT ERSIPHTREI QQALVDAGDK PATFVGSRQW IGSIEVQLVL NQLIGITSKI LFVSQGSEIA SQGRELANHF QSEGTPVMIG GGVL AHTILG VAWNEITGQI KFLILDPHYT GAEDLQVILE KGWCGWKGPD FWNKDAYYNL CLPQRPNMI
Source	Yeast
Target Names	UFSP2
Protein Names	Recommended name: Ufm1-specific protease 2 Short name= UfSP2 EC= 3.4.22.-
Expression Region	1-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.