



Recombinant Human Asparagine synthetase domain-containing protein 1 (ASNSD1)

Product Code	CSB-BP885745HU
Abbreviation	ASNSD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NWL6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	CGICCSVNF SAEHFSQDLK EDLLYNLKQR GPNSSKQLLK SDVNYQCLFS AHVLHLRGVL TTQPVEDERG NVFLWNGEIF SGIKVEAEEN DTQILFNLYS SCKNESEILS LFSEVQGPWS FIYYQASSHY LWFGDRDFGR RLLWHFSNL GKSFCLSSVG TQTSGLANQW QEVPASGLFR IDLKSTVISG CIILQLYPWK YISRENIIEE NVNSLSQISA DLPAFVSVVA NEAKLYLEKP VVPLNMMLPQ AALETHCSNI SNVPPTREIL QVFLTDVHMK EVIQQFIDVL SVAVKKRVLC LPRDENLTAN EVLKTCDRKA NVAILFSGGI DSMVIATLAD RHIPLDEPID LLNVAFIAEE KTMPPTFNRE GNKQKNKCEI PSEEFKSDVA AAAADSPNKH VSVDPDRITGR AGLKELQAVS PSRIWNFVEI NVSMEELQKL RRTRICHLIR PLDVLDDSI GCAVWFASRG IGWLVAQEGV KSYQSNKVV LTGIGADEQL AGYSRHRVRF QSHGLEGLNK EIMMELGRIS SRNLGRDDR IGDHGKEARF PFLDENVVSF LNSLPIWEKA NLTLPRGIGE KLLRLAAVE LGLTASALLP KRAMQFGSRI AKMEKINEKA SDKCGRLQIM SLENLSIEKE TKL
Source	Baculovirus
Target Names	ASNSD1
Protein Names	Recommended name: Asparagine synthetase domain-containing protein 1 Alternative name(s): HCV NS3-transactivated protein 1
Expression Region	2-643
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.