



# Recombinant Human Asparagine synthetase domain-containing protein 1 (ASNSD1)

<b>Product Code</b>	CSB-EP885745HU-B
<b>Abbreviation</b>	ASNSD1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NWL6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	CGICCSVNF SAEHFSQDLK EDLLYNLKQR GPNSSKQLLK SDVNYQCLFS AHVLHLRGVL TTQPVEDERG NVFLWNGEIF SGIKVEAEEN DTQILFNYLS SCKNESEILS LFSEVQGPWS FIYYQASSHY LWFGRDFFGR RLLWHFSNL GKSFCLSSVG TQTSGLANQW QEVPASGLFR IDLKSTVISG CIILQLYPWK YISRENIIEE NVNSLSQISA DLPAFVSVVA NEAKLYLEKP VVPLNMMLPQ AALETHCSNI SNVPPTREIL QVFLTDVHMK EVIQQFIDVL SVAVKKRVLK LPRDENLTAN EVLKTCDRKA NVAILFSGGI DSMVIATLAD RHIPLDEPID LLNVAFIAEE KTMPPTFNRE GNKQKNKCEI PSEEFKSDVA AAAADSPNKH VSVDPDRITGR AGLKELQAVS PSRIWNFVEI NVSMEELQKL RRTRICHLIR PLDVLDDSI GCAVWFASRG IGWLVAQEGV KSYQSNKVV LTGIGADEQL AGYSRHRVRF QSHGLEGLNK EIMMELGRIS SRNLGRDDR IGDHGKEARF PFLDENVVSF LNSLPIWEKA NLTLPRGIGE KLLRLAAVE LGLTASALLP KRAMQFGSRI AKMEKINEKA SDKCGRLQIM SLENLSIEKE TKL
<b>Source</b>	E.coli
<b>Target Names</b>	ASNSD1
<b>Protein Names</b>	Recommended name: Asparagine synthetase domain-containing protein 1 Alternative name(s): HCV NS3-transactivated protein 1
<b>Expression Region</b>	2-643
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.