



# Recombinant Human BRISC and BRCA1-A complex member 1 (BABAM1)

<b>Product Code</b>	CSB-MP885748HU
<b>Abbreviation</b>	BABAM1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NWV8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEVAEPSSPT EEEEEEEHS AEPRPRTRSN PEGAEDRAVG AQASVGSRSE GEGEAASADD GSLNTSGAGP KSWQVPPAP EVQIRTPRVN CPEKVIICLD LSEEMSLPKL ESFNGSKTNA LNVSQMIEM FVRTKHKIDK SHEFALVVVN DDTAWLSGLT SDPRELCSCL YDLETASCST FNLEGLFSLI QKTELPTVE NVQTIIPPYV VRTILVYSRP PCQPQFSLTE PMKKMFQCPY FFDVVYIHN GTEEKEEEMS WKDMFAFMGS LDTKGTSYKY EVALAGPALE LHNCMAKLLA HPLQRPCQSH ASYSLLEEED EAIEVEATV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	BABAM1
<b>Protein Names</b>	Recommended name: BRISC and BRCA1-A complex member 1 Alternative name(s): Mediator of RAP80 interactions and targeting subunit of 40 kDa New component of the BRCA1-A complex
<b>Expression Region</b>	1-329
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.