



Recombinant Human BRISC and BRCA1-A complex member 1 (BABAM1)

Product Code	CSB-YP885748HU
Abbreviation	BABAM1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NWV8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEVAEPSSPT EEEEEEEHS AEPRPRTRSN PEGAEDRAVG AQASVGSRSRSE GEGEAASADD GSLNTSGAGP KSWQVPPAP EVQIRTPRVN CPEKVIICLD LSEEMSLPKL ESFNGSKTNA LNVSQKMIEM FVRTKHKIDK SHEFALVVVN DDTAWLSGLT SDPRELCSCS YDLETASCST FNLEGLFSLI QQKTELPVTE NVQTIPPPYV VRTILVYSRP PCQPQFSLTE PMKKMFQCPY FFFDVVYIHN GTEEKEEEMS WKDMFAFMGS LDTKGTSYKY EVALAGPALE LHNCMAKLLA HPLQRPCQSH ASYSLLEEED EAIEVEATV
Source	Yeast
Target Names	BABAM1
Protein Names	Recommended name: BRISC and BRCA1-A complex member 1 Alternative name(s): Mediator of RAP80 interactions and targeting subunit of 40 kDa New component of the BRCA1-A complex
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.