



Recombinant Human Mediator of RNA polymerase II transcription subunit 29 (MED29)

Product Code	CSB-MP885752HU
Abbreviation	MED29
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NX70
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAASQQQASA ASSAAGVSGP SSAGGPGPQQ QPQPPAQLVG PAQSGLLQQQ QQDFDPVQRY KMLIPQLKES LQTLMKVAAQ NLIQNTNIDN GQKSSDGPIQ RFDKCLEEFY ALCDQLELCL RLAHECLSQS CDSAKHSPTL VPTATKPDV QPDSLPPYQY LAVIKAQISC AKDIHTALLD CANKVTGKTP APPAGPGGTL
Source	Mammalian cell
Target Names	MED29
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 29 Alternative name(s): Intersex-like protein Mediator complex subunit 29
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.