



Recombinant Human Probable alpha-ketoglutarate-dependent dioxygenase ABH4 (ALKBH4)

Product Code	CSB-BP885759HU
Abbreviation	ALKBH4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NXW9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AAAAAETPE VLRECGCKGI RTCLICERQR GSDPPWELPP AKTYRFIYCS DTGWAVGTEE SDFEGWAFPF PGVMLIEDFV TREEEAELVR LMDRDPWKLS QSGRRKQDYG PKVNFRRKQKL KTEGFCGLPS FSREVVRRMG LYPGLEGFRP VEQCNDYCP ERGSAIDPHL DDAWLWGERL VSLNLLSPTV LSMCREAPGS LLLCSAPSAA PEALVDSVIA PSRSVLCQEV EVAIPLPARS LLVLTGAARH QWKHAIHRRH IEARRVCVTF RELSAEFGPG GRQQELGQEL LRIALSFQGR PV
Source	Baculovirus
Target Names	ALKBH4
Protein Names	Recommended name: Probable alpha-ketoglutarate-dependent dioxygenase ABH4 EC= 1.14.11.- Alternative name(s): Alkylated DNA repair protein alkB homolog 4
Expression Region	2-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.