



Recombinant Human Core histone macro-H2A.2 (MACROH2A2)

Product Code	CSB-YP885785HU
Abbreviation	H2AFY2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P0M6
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSGRSKGKKM SKLSRSARAG VIFPVGRLMR YLKKGTFKYR ISVGAPVYMA AVIEYLAAEI LELAGNAARD NKKARIAPRH ILLAVANDEE LNQLLKGVTI ASGGVLPRIH PELLAKKRGT KGKSETILSP PPEKRGRKAT SGKKGGKSK AAKPRTSKKS KPKDSDKEGT SNSTSEDGPG DGFTILSSKS LVLGQKLSLT QSDISHIGSM RVEGIVHPTT AEIDLKEDIG KALEKAGGKE FLETVKELRK SQGPLEVAEA AVSQSSGLAA KFVIHCHIPQ WGSDKCEEQL EETIKNCLSA AEDKCLKSVA FPPFPSGRNC FPKQTAAQVT LKAISAHFDD SSASSLKNVY FLLFDESEIG IYVQEMAKLD AK
Source	Yeast
Target Names	H2AFY2
Protein Names	Recommended name: Core histone macro-H2A.2 Short name= Histone macroH2A2 Short name= mH2A2
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.