



Recombinant Human Protein RCC2 (RCC2)

Product Code	CSB-BP885796HU
Abbreviation	RCC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P258
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPRKKAAAAA WEEPSSGNGT ARAGPRKRGG PAGRKREPERE RCSSSSGGGS SGDEDGLELD GAPGGGKRAA RPATAGKAGG AAVVITEPEH TKERVKLEGS KCKGQLLIFG ATNWDLIGRK EVPKQQAAYR NLGQNLWGPH RYGCLAGVRV RTVVS GSCAA HSLITTEGK LWSWGRNEKG QLGHGDTKRV EAPRLIEGLS HEVIVSAACG RNHTLALTET GSVFAFGENK MGQLGLGNQT DAVPSPAQIM YNGQPITKMA CGAEFSMIMD CKGNLYSFGC PEYGQLGHNS DGKFIARAQR IEYDCELVPR RVAIFIEKTK DGQILVPPNV VVRDVACGAN HTLVLDSQKR VFSWGFGGYG RLGHAEQKDE MVPRLVKLFD FPGRGASQIY AGYTCSFAVS EVGGLFFWGA TNSRESTMY PKAVQDLGW RIRSLACGKS SIIVADEST ISWGPSPTFG ELGYGDHKPK SSTA AQEVKT LDGIFSEQVA MGYSHSLVIA RDESETEKEK IKKLPEYNPR TL
Source	Baculovirus
Target Names	RCC2
Protein Names	Recommended name: Protein RCC2 Alternative name(s): RCC1-like protein TD-60 Telophase disk protein of 60 kDa
Expression Region	1-522
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.