



# Recombinant Human Protein RCC2 (RCC2)

<b>Product Code</b>	CSB-EP885796HU-B
<b>Abbreviation</b>	RCC2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9P258
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPRKKA AAAA WEEPSSGNGT ARAGPRKRGG PAGRKREPE RCSSSSGGGS SGDEDGLELD GAPGGGKRAA RPATAGKAGG AAVVITEPEH TKERVKLEGS KCKGQLLIFG ATNWDLIGRK EVPKQQAAYR NLGQNLWGPH RYGCLAGVRV RTVVS GSCAA HLLITTEGK LWSWGRNEKG QLGHGDTKRV EAPRLIEGLS HEVIVSAACG RNHTLALTET GSVFAFGENK MGQLGLGNQT DAVPSPAQIM YNGQPITKMA CGAEFSMIMD CKGNLYSFGC PEYGQLGHNS DGKFIARAQR IEYDCELVPR RVAIFIEKTK DGQILVPPNV VVRDVACGAN HTLVLDSQKR VFSWGFGGYG RLGHAEQKDE MVPRLVKLFD FPGRGASQIY AGYTCSFAVS EVGGLFFWGA TNSRESTMY PKAVQDLGW RIRSLACGKS SIIVAADEST ISWGPSPTFG ELGYGDHKPK SSTA AQEVKT LDGIFSEQVA MGYSHSLVIA RDESETEKEK IKKLPEYNPR TL
<b>Source</b>	E.coli
<b>Target Names</b>	RCC2
<b>Protein Names</b>	Recommended name: Protein RCC2 Alternative name(s): RCC1-like protein TD-60 Telophase disk protein of 60 kDa
<b>Expression Region</b>	1-522
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.