



Recombinant Human Intraflagellar transport protein 80 homolog (IFT80)

Product Code	CSB-BP885802HU
Abbreviation	IFT80
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P2H3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRLKISLLKE PKHQELVSCV GWTTAEELYS CSDDHQIVKW NLLTSETTQI VKLPDDIYPI DFHWFPKSLG VKKQTQAESF VLTSSDGKFH LISKLGRVEK SVEAHCGAVL AGRWNYEGTA LVTVGEDGQI KIWSKTGMLR STLAQQGTPV YSVAWGPDSE KVLTYTAGKQL IIKPLQPNAK VLQWKAHDGI ILKVDWNSVN DLILSAGEDC KYKVWDSYGR PLYNSQPHEH PITSAWAPD GELFAVGSFH TLRLCDKTGW SYALEKPNTG SIFNIAWSID GTQIAGACGN GHVVFHVVVE QHWEWKNFQV TLTKRRAQMV RNVLNDAVDL LEFRDRVIKA SLNYAHLVVS TSLQCYVFST KNWNTPIIFD LKEGTVSLIL QAERHLLVD GSSIYLYSYE GRFISSPKFP GMRTDILNAQ TVSLSNDTIA IRDKADEKII FLFEASTGKP LGDGKFLSHK NEILEIALDQ KGLTNDRKIA FIDKNRDL CI TSVKRFSGKEE QIIKLGTMVH TLAWNDCNI LCGLQDTRFI VWYYPNTVYV DRDILPKTLY ERDASEFSKN PHIVSFVGNQ VTIRRADGSL VHSITPYPA ILHEYVSSSK WEDAVRLCRF VKEQTMWACL AAMAVANRDM TTAEIAYAAI GEIDKVQYIN SIKNLPSKES KMAHILLFSG NIQEAIEVLL QAGLVYQAIQ ININLYNWER ALELAVKYKT HVDTVLAYRQ KFLETFGKQE TNKRYLHYAE GLQIDWEKIK AKIEMEITKE REQSSSSQSS KSIGLKP
Source	Baculovirus
Target Names	IFT80
Protein Names	Recommended name: Intraflagellar transport protein 80 homolog Alternative name(s): WD repeat-containing protein 56
Expression Region	1-777
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.