



Recombinant Danio rerio Protein L-Myc-1a (mycl1a)

Product Code	CSB-EP886088DIL-B
Abbreviation	mycl1a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9PSI9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MPGVNTQTLY GWDMEHYFYD EMDTEEDFFK STAPSEDIWK KFELLPTPPM SPSRTLGDGW LFPIPGDRLG WAPPKVLTCDEEYEGHLHKFD PLDIFGNLGS IVIKDCMWSG LSTSHRLEKV AHGERAPVAA LIQNSTAQKA ARAASGTPVT GSQAAQCVSP AAVLELPVPH KKVAAGSSGS ECRSDSSDDD EDDDEIDVVT VDNRPKRGRP PSRRTPVTIT VSADPFGPCP KRFHVSLHRQ QHNYAAPSPD TDPEDDFEIE PVSKRPRLES SSAPSSPLSS PATSDSEDST EQRRNFLERK RRDDLRSRFQ ALREEIPGLS GSSKTSKVAI LTQATDYLLQ LHSSQRRQAQ EKRKLKAKQQ QLLRRISALQ NS
Source	E.coli
Target Names	mycl1a
Protein Names	Recommended name: Protein L-Myc-1a Alternative name(s): Protein L-Myc 1 Short name= zL-Myc
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.