



Recombinant Danio rerio Protein NipSnap homolog 2 (gbas)

Product Code	CSB-MP886100DIL
Abbreviation	gbas
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9PU58
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MATRVLHSSC SGLYRAAGPA RGKGHATAVI RSLSASHNRP REDSWFKSLF VRKVDPRKDA HSHLLAKKED NNLYKIQFHN VKPECLEAYN KLCEDVLTNI HTDKAYPCEL VGTWNTWYGE QDQAVHLWRY RGGYPALTEV MSKLNKKEF LEYRSEKGM LLSRRNQLLL EFSFWNEPVP RDGPNIYELR SYQLRPGTMI EWGNYWARAI GYRQHNREAV GGFFSQIGDL YMVHHLWAYK DLQSREDTRN AAWQHEGWDE VVYYTVPLIQ HMESRIMIPL KNSPLK
Source	Mammalian cell
Target Names	nipsnap2
Protein Names	Recommended name: Protein NipSnap homolog 2 Short name= NipSnap2 Alternative name(s): Glioblastoma-amplified sequence
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.