



Recombinant Human immunodeficiency virus type 1 group M subtype K Gag-Pol polyprotein (gag-pol), partial

| | |
|--------------------------|---|
| Product Code | CSB-BP886136HLQ |
| Abbreviation | gag-pol |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9QBZ9 |
| Product Type | Recombinant Protein |
| Immunogen Species | Human immunodeficiency virus type 1 group M subtype K (isolate 97ZR-EQTB11) (HIV-1) |
| Purity | ≥85% (SDS-PAGE) |
| Source | Baculovirus |
| Target Names | gag-pol |
| Protein Names | Recommended name: Gag-Pol polyprotein Alternative name(s): Pr160Gag-Pol Cleaved into the following 11 chains: 1. Matrix protein p17 Short name= 2. MA 3. Capsid protein p24 Short name= 4. CA 5. Spacer peptide p2 6. Nuc |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |