



Recombinant Arabidopsis thaliana Probable glutamate dehydrogenase 3 (GSH3)

Product Code	CSB-EP886609DOA-B
Abbreviation	GSH3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9S7A0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNALAATNRN FKLASRLGL DSKLEKSLLI PFREIKVECT IPKDDGTLAS FVGFRVQHDN ARGPMKGGIR YHPEVEPDEV NALAQLMTWK TAVAKIPYGG AKGGIGCDPS ELSLSELERL TRVFTQKIHD LIGIHTDVPA PDMGTGPQTM AWILDEYSKF HGHSPAVVTG KPIDLGGSLG RDAATGRGVL FATEALLNEH GKTISGQRFA IQGFGNVGSW AAKLISDKGG KIVAVSDVTG AIKNNNGIDI LSLLEHAEEN RGIKGF DGAD SIDPDSILVE DCDILVPAAL GGVINRENAN EIKAKFII EG ANHPTDPEAD EILKKKG VMI LPDIYANSGG VTVSYFEWVQ NIQGF MWDEE KVNRELKTYM TRGFKDLKEM CQTHSCDLRM GAFTLGINRV AQATTIRGWG S
Source	E.coli
Target Names	GSH3
Protein Names	Recommended name: Probable glutamate dehydrogenase 3 Short name= GDH 3 EC= 1.4.1.3
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.