



Recombinant Arabidopsis thaliana Putative DNA-binding protein ESCAROLA (ESC)

Product Code	CSB-MP886610DOA
Abbreviation	ESC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9S7C9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEGGYEQGGG ASRYFHNLFR PEIHHQQLQP QGGINLIDQH HHQHQQHQQQ QQPSDDSRES DHSNKDHHQQ GRPDSDPNTS SSAPGKRPRG RPPGSKNKAK PPIIVTRDSP NALRSHVLEV SPGADIVESV STYARRRGRG VSVLGGNGTV SNVTLRQPVT PGNGGGVSGG GGVTLHGRF EILSLTGTVL PPPAPPGAGG LSIFLAGGGQ QVVGGSVVAP LIASAPVILM AASFNAVFE RLPIEEEEEE GGGGGGGGGG GPPQMQQAPS ASPPSGVTGQ GQLGGNVGGY GFSGDPHLLG WGAGTPSRPP F
Source	Mammalian cell
Target Names	AHL27
Protein Names	Recommended name: Putative DNA-binding protein ESCAROLA
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.