



# Recombinant Arabidopsis thaliana Nitrile-specifier protein 1 (NSP1)

<b>Product Code</b>	CSB-BP886656DOA
<b>Abbreviation</b>	NSP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9SDM9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAQKLEAKGG EMGDVWDDGV YENVRKVYVG QAQYGIAFVK FEYVNGSQVV VGDEHGKKTE LGVEEFEIDA DDYIVYVEGY REKVNDMTSE MITFLSIKTF KGKTSHPTEK RPGVKFVLHG GKIVGFHGRS TDVLHSLGAY VLSSTIKLL GKWIKVEQKG EGPGLRCSHG IAQVGNKIYS FGGEFTPNQP IDKHLYVFDL ETRTWSISPA TGDVPHLSCL GVRMVSVGST LYVFGGRDAS RQYNGFYSFD TTTNEWKLLT PVEEGTPRS FHSMAADEEN VYVFGGVSAT ARLNTLDSYN IVDKKWFHCS TPGDSLARG GAGLEVQVK VVVYGFNGC EVDDVHYYDP VQDKWTQVET FGVRPERSV FASAAIGKHI VIFGGEIAMD PLAHVGPQL TDGTFALDTE TLQWERLDKF GGEETPSSR GWTASTTATI DGKKGLVMHG GKAPTNDKFD DLFFYGIDSA
<b>Source</b>	Baculovirus
<b>Target Names</b>	NSP1
<b>Protein Names</b>	Recommended name: Nitrile-specifier protein 1 Short name= AtNSP1
<b>Expression Region</b>	1-470
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.