



Recombinant Arabidopsis thaliana Nitrile-specifier protein 1 (NSP1)

Product Code	CSB-MP886656DOA
Abbreviation	NSP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SDM9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAQKLEAKGG EMGDVWDDGV YENVRKVYVG QAQYGIAFVK FEYVNGSQVV VGDEHGKKTE LGVEEFEIDA DDYIVYVEGY REKVNDMTSE MITFLSIKTF KGKTSHPTEK RPGVKFVLHG GKIVGFHGRS TDVLHSLGAY VLSSTIKLL GKWIKVEQKG EGPGLRCSHG IAQVGNKIYS FGGEFTPNQP IDKHLYVFDL ETRTWSISPA TGDVPHLSCL GVRMVS VGST LYVFGGRDAS RQYNGFYSFD TTTNEWKLLT PVEEGTPRS FHSMAADEEN VYVFGGVSAT ARLNTLDSYN IVDKKWFHCS TPGDSL TARG GAGLEV VQ GK VVVYGFNGC EVDDVHYYDP VQDKWTQVET FGV RPSERSV FASAAIGKHI VIFGGEIAMD PLAHVGPQQL TDGTFALDTE TLQWERLDFK G GEEETPSSR GWTASTTATI DGKKGLVMHG GKAPTND RFD DLFFYGIDSA
Source	Mammalian cell
Target Names	NSP1
Protein Names	Recommended name: Nitrile-specifier protein 1 Short name= AtNSP1
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.