



Recombinant Arabidopsis thaliana MADS-box transcription factor ANR1 (ANR1)

Product Code	CSB-MP886685DOA
Abbreviation	ANR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SI38
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGRGKIVIRR IDNSTSRQVT FSKRRSGLLK KAKELSILCD AEVGVIIFSS TGKLYDYASN SSMKTIERY NRVKEEQHQL LNHASEIKFW QREVASLQQQ LQYLQECHRK LVGEELSGMN ANDLQNLEDQ LVTSLKGVRL KKDQLMTNEI RELNRKGQII QKENHELQNI VDIMRKENIK LQKKVHGRTN AIEGNSSVDP ISNGTTTYAP PQLQLIQLQP APREKSIRLG LQLS
Source	Mammalian cell
Target Names	ANR1
Protein Names	Recommended name: MADS-box transcription factor ANR1 Alternative name(s): Protein AGAMOUS-LIKE 44 Protein ARABIDOPSIS NITRATE REGULATED 1
Expression Region	1-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.