



Recombinant Arabidopsis thaliana Probable WRKY transcription factor 17 (WRKY17)

Product Code	CSB-YP886700DOA
Abbreviation	WRKY17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9SJA8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTVDIMRLPK MEDQTAIQEA ASQGLKSMEH LIRVLSNRPE ERNVDCSEIT DFTVSKFKKV ISLLNRSGHA RFRRGPVHSP PSSSVPPPVK VTPAPTQIS APAPVSFVQA NQQSVTLDFR RPSVFGAKTK SSEVVEFAKE SFSVSSNSSF MSSAITGDGS VSKGSSIFLA PPAVPVTSS GKPLSGLPY RKRCFEHDHS EGFSGKISGS GNGKCHCKKS RKNRMKRTVR VPAVSAKIAD IPPDEYSWRK YGQKPIKGSP HPRGYKYST FRGCPARKHV ERALDDSTML IVTYEGEHRH HQSTMQEHTV PSVSGLVFGS A
Source	Yeast
Target Names	WRKY17
Protein Names	Recommended name: Probable WRKY transcription factor 17 Alternative name(s): WRKY DNA-binding protein 17
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.